

MACSQuant® Flow Cytometers

Engineering the future of automation





At the heart of intelligent design

The first MACSQuant® Analyzer was engineered to enable demanding and diverse flow cytometry applications of researchers around the world.

At Miltenyi Biotec, we have since continued to expand our instrument portfolio to meet the evergrowing complexities in cell analysis. Innovation is at the core of each MACSQuant Instrument, offering automated flow cytometry solutions for a diverse array of appliations.

Wondering which MACSQuant Flow Cytometer is the best fit for your lab? Explore your options in the MACSQuant family of flow cytometers.

- MACSQuant Analyzer 10
 automating and standardizing key assays
- MACSQuant Analyzer 16
 increasing data acquisition parameters with
 more colors
- MACSQuant X
 high-throughput applications from
 streamlined fluidics and sample processing
- MACSQuant VYB
 red fluorescent protein analysis and
 optimized FITC and PE applications

	•.	.	0.	<u>.</u>
	• MACS 10	MACS Quant	MACS 16	MACS QuantVYB
Automation features	•	• Quality (•	•
Automated acquisition	•	•	•	•
Automated labeling	•	•	•	•
Automated start up, shutdown, cleaning	•	•	•	•
Automated compensation	•	•	•	•
Optics				
10 parameter	•	•		•
16 parameter			•	
405 nm laser	•	•	•	•
488 nm laser	•	•	•	•
561 nm laser				•
640 nm laser	•	•	•	
Violet 405 nm side scatter			•	
Sample processing				
Needle-arm vibration mixing		•		
Orbital sample mixing		•		
Pipette sample mixing	•		•	•
24-tube rack processing	•	•	•	•
96-well plate processing	•	•	•	•
384-well plate processing		•		
High-throughput plate screening		•		
Sample cooling (chill racks)	•		•	•
Low carryover: standard mode < 0.01%	•	•	•	•
Unique features				
Absolute cell counting	•	•	•	•
Rare cell enrichment	•		•	•
Remote support	•	•	•	•
Multi-instrument alignment	•	•	•	•
21 CFR Part 11 compliant software	•	•	•	•





MACSQuant® Analyzer 10

Your trusted partner in automated flow cytometry

The MACSQuant® Analyzer 10 brings to life the versatility and power required for modern flow cytometry applications. Whether mining for rare cells, analyzing the efficiency of your cell manufacturing process or investigating signaling pathways, you are equipped for the task at hand.


Truly hands-free operation

With a range of automated features, the MACSQuant Analyzer 10 lays the foundation for true automation present in all instruments belonging to the MACSQuant Family of flow cytometers.



Application highlight: Automated assays

- In-process and quality control
- Cellular therapy, i.e., chimeric antigen receptor (CAR) T cell manufacturing
- Toxicological assays
- · Whole blood analysis
- Cell counting
- Rare cell analysis
- Small particles/exosomes
- Cytokine arrays



MACSQuant Analyzer 10 optical bench

405 nm

450/50 nm VioBlue® 525/50 nm VioGreen™ 488 nm

488/10 nm	FSC
488/10 nm	SSC
525/50 nm	FITC
585/40 nm	PE
655-730 nm	PerCP-Vio® 700
750 nm LP	PE-Vio® 770

655–730 nm	APC
750 nm l P	APC-Vio® 770




MACSQuant® X


Reliable high-throughput flow cytometry


Engineered from the ground up as a compact flow cytometer, the MACSQuant® X was made to provide speed and reliability in high-throughput screening. It allows you to get more data, in less time, with less effort. Perfect for laboratories that want to maximize efficiency, while enjoying the certainty of getting reliable results from every tube, 96- or 384-well plate. The MACSQuant X has been designed to meet high-throughput cell analysis needs for the most challenging large-scale screening settings.

Technology highlight: Integration with liquid handling systems

- Robust flow cytometer capable of processing up to 9,200 samples per day (96-well plate within 15 minutes, 384-well plate in under 60 minutes)
- Seamless integration into a variety of liquid handling systems







MACSQuant X optical bench

405 nm

450/50 nm	VioBlue®
525/50 nm	VioGreen™

488 nm

488/10 nm	FSC
488/10 nm	SSC
525/50 nm	FITC
585/40 nm	PE
655–730 nm	PerCP-Vio® 700
750 nm LP	PE-Vio® 770

655–730 nm	APC
750 nm LP	APC-Vio® 770



MACSQuant® Analyzer 16

Colorful new possibilities in automated flow cytometry

Expanding fully automated data acquisition to 16-parameters is what the MACSQuant® Analyzer 16 is all about. The instrument's compact design is suited for basic research as well as advanced immune monitoring applications. It provides users with the flexibility and customization known from the MACSQuant Analyzer 10, with an additional six fluorescent channels to meet the increasing demands of modern laboratories.

Unique features

- Expanded fluorescent channels: Amplify your data with up to 16 parameters, while using 40% less sample material
- Violet side scatter: Enhance your small particle analysis with the optimal laser line.

Technology highlight: More acquisition parameters and multi-instrument alignment

- Explore more cell subsets and their activation status by utilizing a broader range of fluorochrome and antibody options
- Improved consistency and reproducibility by easy aligning of instrument settings for intra- and inter-lab assays or applications.
 Down the hall or around the world!



MACSQuant Analyzer 16 optical bench

405 nm

405/10 nm	SSC
450/50 nm	VioBlue®
525/50 nm	VioGreen™
579/34 nm	BV 570™
615/20 nm	BV 605™
667/30 nm	BV 650™

488 nm

488/10 nm	FSC
525/50 nm	FITC
579/34 nm	PE
615/20 nm	PE-Vio® 615
667/30 nm	PerCP
725/40 nm	PerCP-Vio® 700
785/62 nm	PE-Vio® 770

667/30 nm	APC
725/40 nm	Alexa Fluor® 700
785/62 nm	APC-Vio® 770





MACSQuant® VYB

Your benchmark for fluorescent protein flow cytometry

On top of the performance, automation, compactness and convenience associated with the MACSQuant® Family of flow cytometers, the MACSQuant VYB comes with a uniquely configured optical bench, featuring violet, yellow, and blue lasers. It is a perfect match for labs utilizing fluorescent protein reporters or for labs that want to simply use PE and FITC conjugates.

Technology highlight: Finding the right fit for fluorescent proteins with the yellow laser

- Fully utilize fluorescent proteins and reporter applications that require the Violet, Blue or Yellow laser such as mCherry, GFP and CFP (figure 2)
- Optimal excitation of PE and PE tandem dyes
- Excite and analyze APC and APC tandems
- Minimize compensation for FITC/GFP and PE conjugates (figure 3)




Figure 2: The MACSQuant VYB can detect a wide range of fluorescent proteins. Histograms of cells transfected with CFP, GFP, YFP (top row), mCherry, and mKate 2 (bottom row) are easily distinguished from nontransfected cells (shown in purple). GFP and YFP are often difficult to distinguish using flow cytometry due to similarities in their emission spectra. With the MACSQuant VYB, YFP (shown in yellow) and GFP (shown in green) cells are clearly distinguishable when observed in the PI channel (bottom right).




Figure 3: Cells expressing GFP were stained using a PE-conjugated antibody against a dim marker and analyzed using the MACSQuant VYB. Excitation of both PE and GFP with the blue 488 nm laser only minimally resolves the two signals (**A**). In comparison, excitation of PE by the yellow 561 nm laser and of GFP by the blue 488 nm laser fully resolves these two signals without the need for any compensations (**B**).

MACSQuant VYB optical bench

405 nm

450/50 nm	VioBlue®, BFP
525/50 nm	VioGreen™

488 nm

525/50 nm	GFP, FITC
614/50 nm	PI, LSS-mKate

561/10 nm	FSC
561/10 nm	SSC
586/15 nm	PE
615/20 nm	mCherry, PE-Vio® 615
655–730 nm	mKate, APC
750 nm LP	PE-Vio® 770, APC-Vio® 770

▶ miltenyibiotec.com/macsquant



Germany/Austria/

Switzerland Miltenyi Biotec GmbH

Friedrich-Ebert-Straße 68 51429 Bergisch Gladbach Germany Phone +49 2204 8306-0 Fax +49 2204 85197 macs@miltenyibiotec.de

USA/Canada

Miltenyi Biotec Inc. 2303 Lindbergh Street Auburn, CA 95602, USA Phone 800 FOR MACS Phone +1 530 888 8871 Fax +1 877 591 1060 macs@miltenyibiotec.com

Miltenyi Biotec Australia Pty. Ltd. Unit 16A, 2 Éden Park Drive Macquarie Park NSW 2113 Australia Phone +61 2 8877 7400 Fax +61 2 9889 5044 macs@miltenyibiotec.com.au

Benelux

Miltenyi Biotec B.V. Schipholweg 68 H 2316 XE Leiden The Netherlands macs@miltenyibiotec.nl Customer service

The Netherlands

Phone 0800 4020120 Fax 0800 4020100

Customer service Belgium Phone 0800 94016 Fax 0800 99626

Customer service Luxembourg Phone 800 24971 Fax 800 24984

China

Miltenyi Biotec Technology & Trading (Shanghai) Co., Ltd. Rooms 2303 and 2309 No. 319, Xianxia Road **Changning District** 200051 Shanghai, P.R. China Phone +86 21 62351005 Fax +86 21 62350953 macs@miltenyibiotec.com.cn

Miltenyi Biotec SAS 10 rue Mercoeur 75011 Paris, France Phone +33 1 56 98 16 16 Fax +33 1 56 98 16 17 macs@miltenyibiotec.fr

Miltenyi Biotec S.r.l. Via Paolo Nanni Costa, 30 40133 Bologna Italy Phone +39 051 6 460 411 Fax +39 051 6 460 499 macs@miltenyibiotec.it

Japan

Miltenyi Biotec K.K. Nittsu-Eitai Building 5F 16-10 Fuyuki, Koto-ku, Tokyo 135-0041, Japan Phone +81 3 5646 8910 Fax +81 3 5646 8911 macs@miltenyibiotec.jp

Nordics and Baltics

Miltenyi Biotec Norden AB Scheelevägen 17 223 70 Lund Sweden macs@miltenyibiotec.se Customer service Sweden

Phone 0200-111 800 Fax 046-280 72 99

Customer service Denmark Phone 80 20 30 10 Fax +46 46 280 72 99

Customer service Norway, Finland, Iceland, and Baltic countries Phone +46 46 280 72 80 Fax +46 46 280 72 99

Singapore

Miltenyi Biotec Asia Pacific Pte Ltd. 100 Beach Road #28-06 to 28-08 Shaw Tower Singapore 189702 Phone +65 6238 8183 Fax +65 6238 0302 macs@miltenyibiotec.com.sg

South Korea

Miltenyi Biotec Korea Co., Ltd Arigi Bldg. 8F 562 Nonhyeon-ro Gangnam-gu Seoul 06136, South Korea Phone +82 2 555 1988 Fax +82 2 555 8890 macs@miltenyibiotec.co.kr

Spain

Miltenyi Biotec S.L. C/Luis Buñuel 2 Ciudad de la Imagen 28223 Pozuelo de Alarcón (Madrid) Phone +34 91 512 12 90 Fax +34 91 512 12 91 macs@miltenyibiotec.es

United Kingdom

Miltenyi Biotec Ltd. Almac House, Church Lane Bisley, Surrey GU24 9DR, UK Phone +44 1483 799 800 Fax +44 1483 799 811 macs@miltenyibiotec.co.uk

www.miltenyibiotec.com

Miltenyi Biotec provides products and services worldwide. Visit www.miltenyibiotec.com/local to find your nearest Miltenyi Biotec contact.